

# Chapter 12

## Identification of Open Chromatin Regions in Plant Genomes Using ATAC-Seq

Marko Bajic, Kelsey A. Maher, and Roger B. Deal

### Abstract

Identifying and characterizing highly accessible chromatin regions assists in determining the location of genomic regulatory elements and understanding transcriptional regulation. In this chapter, we describe an approach to map accessible chromatin features in plants using the Assay for Transposase-Accessible Chromatin, combined with high-throughput sequencing (ATAC-seq), which was originally developed for cultured animal cells. This technique utilizes a hyperactive Tn5 transposase to cause DNA cleavage and simultaneous insertion of sequencing adapters into open chromatin regions of the input nuclei. The application of ATAC-seq to plant tissue has been challenging due to the difficulty of isolating nuclei sufficiently free of interfering organellar DNA. Here we present two different approaches to purify plant nuclei for ATAC-seq: the INTACT method (Isolation of Nuclei TAgged in specific Cell Types) to isolate nuclei from individual cell types of the plant, and tissue lysis followed by sucrose sedimentation to isolate sufficiently pure total nuclei. We provide detailed instructions for transposase treatment of nuclei isolated using either approach, as well as subsequent preparation of ATAC-seq libraries. Sequencing-ready ATAC-seq libraries can be prepared from plant tissue in as little as one day. The procedures described here are optimized for *Arabidopsis thaliana* but can also be applied to other plant species.

**Key words** ATAC-seq, INTACT system, Chromatin, Nucleus, Transposition, Nucleosome, Transcription factor, Enhancer

---


### 1 Introduction

Plants are sessile organisms that must precisely regulate their transcription in response to their environment, as well as for proper development, growth, and homeostasis. Transcription is associated with regions of relatively open chromatin, in which *cis*-regulatory elements such as enhancers and promoters can recruit transcription factors and RNA polymerase II to transcribe DNA [1]. Binding of transcription factors to DNA generally results in the depletion of nucleosomes, rendering these regions hypersensitive to nucleases. Characterizing such regulatory regions throughout the genome has therefore relied on methods that combine enzymatic digestion of nuclear DNA and high-throughput sequencing, such as

micrococcal nuclease sequencing (MNase-seq, *see* Chapter 10 [2]) and DNase I hypersensitivity sequencing (DNase-seq) [3, 4]. Alternatively, regulatory regions can be inferred by chromatin immunoprecipitation sequencing (ChIP-seq, *see* Chapter 5 [5]) where antibodies are used to pull down transcription factors or histone marks associated with active transcription [6].

An improved method for identifying accessible regions of chromatin and transcription factor binding is the Assay for Transposase-Accessible Chromatin with high-throughput sequencing (ATAC-seq) [7, 8]. This method uses a hyperactive Tn5 transposase to integrate preloaded sequencing adapters into regions of open chromatin (Fig. 1a). ATAC-seq is a fast protocol with simple library amplification steps and requires very small amounts of starting material, making it a vast improvement over alternative methods. However, a drawback of this protocol is that the hyperactive Tn5 transposase also targets sources of extranuclear genetic material, including the genomes of mitochondria and chloroplasts. This decreases the proportion of reads that map to the nuclear genome, reducing the amount of information that can be used to identify regulatory regions of open chromatin. Such extranuclear reads must be discarded at the start of the data analysis process, diminishing the efficiency of the assay both in terms of cost and in effective use of materials. To gain the maximum efficiency of this powerful procedure, input material free from extranuclear genetic material, such as purified nuclei, is the ideal input for ATAC-seq.

In this chapter, we describe the use of two different methods to isolate either total nuclei from tissues or nuclei from specific cell types of *Arabidopsis thaliana* (Fig. 1b). To isolate total nuclei from plant tissue we use extraction buffers with a non-ionic detergent to lyse organelles, followed by sucrose sedimentation to further purify the nuclei [9]. This method of nuclei isolation can be done in any lab on most plant tissues. However, these partially purified nuclei still contain some organellar DNA in addition to nuclear DNA, which reduces the efficiency of Tn5 transposition and results in fewer sequencing reads that map to nuclear DNA. In addition, we describe the Isolation of Nuclei Tagged in specific Cell Types (INTACT) method to isolate nuclei from tissue or from specific cell types [10]. This system uses two transgenes for nuclear targeting for affinity purification: (1) the Nuclear Tagging Fusion (NTF) construct, which encodes a fusion of WPP nuclear envelope-targeting domain, a green fluorescent protein (GFP), and the biotin ligase recognition peptide (BLRP); and (2) an *E. coli* biotin ligase (BirA), which biotinylates the BLRP tag. The BirA is expressed from a constitutive promoter while the NTF is expressed either from a constitutive or cell type-specific promoter. The specificity of the NTF promoter determines which cell types will have biotinylated nuclei, and can then be isolated by affinity purification with streptavidin-coated magnetic beads [11]. A key advantage of the



**Fig. 1** ATAC-seq profiling using nuclei isolated by INTACT or sucrose sedimentation. **(a)** Overview of the ATAC-seq procedure. Nuclei are incubated with sequencing adapter-loaded Tn5 transposase, which diffuses into the nucleus to interact with chromatin. Sequencing adapters are inserted into open chromatin regions, and the fragmented DNA is amplified wherever the sequencing adapters were inserted. This generates a library of DNA fragments in which each end represents an insertion site. The amplified libraries are purified and sequenced with next-generation sequencing. **(b)** Two different methods for purifying nuclei from *Arabidopsis* can be used: (1) INTACT for isolating nuclei from specific cell types, and (2) sucrose sedimentation to isolate total nuclei from input tissue. The two methods have the same initial steps: tissue is collected from a specific part of the plant (root, leaf, or the entire plant), ground to a fine powder, resuspended, filtered, and centrifuged to pellet nuclei and cellular debris. Nuclei isolation using tissue that expresses INTACT transgenes uses streptavidin coated magnetic beads to affinity-purify biotinylated nuclei out of the resuspended pellet. This allows for the isolation of nuclei from specific cell types that express the Nuclear Tagging Fusion (NTF) protein and the biotin ligase BirA, resulting in very low contamination by organellar genomes. Alternatively, total nuclei can be isolated from tissue by resuspending the nuclei/debris pellet in a buffer with Triton X-100 to lyse organelles and centrifuging through a dense sucrose layer. Nuclei isolated from both procedures are stained with DAPI and quantified using a hemocytometer. **(c)** Fluorescent microscope images of nuclei (white arrows) stained with the DNA-binding dye DAPI (blue) isolated either through INTACT or sucrose sedimentation. INTACT isolated nuclei are identified by their DAPI fluorescence and binding to multiple beads (white arrowhead). Beads are easily visualized by increasing the transmission of white light while viewing the nuclei in the DAPI channel. Sucrose sedimentation isolated nuclei (white arrows) are DAPI-stained objects around 4–6  $\mu\text{m}$  in diameter, although they can vary in size and shape depending on starting tissue. Much more cellular debris (white asterisk) is observed in sucrose sedimentation-isolated nuclei as compared to INTACT-purified nuclei, but this should not impact the procedure described here. Each picture contains a 50  $\mu\text{m}$  scale bar shown at the bottom

INTACT approach is not only that the isolated nuclei have less organellar DNA contamination, but also that this method can be used to selectively isolate nuclei from specific cell types. While INTACT is a powerful technique, it does require that stable transgenic lines containing BirA and NTF cassettes for the cell type of interest are available, which are time consuming to generate and can be limiting for many species. Even so, the protocol described here, particularly ATAC-seq using sucrose sedimentation-purified nuclei, can readily be adapted for chromatin profiling in any plant species.

---

## 2 Materials

### 2.1 Equipment

1. Porcelain 50 mL mortar and pestle, or equivalent.
2. Liquid nitrogen.
3. Metal lab spoon.
4. Magnetic rack for 1.5 mL tubes (e.g., DynaMag 2, Life Technologies).
5. Magnetic rack for 15 mL tubes (e.g., DynaMag 15, Life Technologies).
6. 96-well magnetic separator plate (e.g., MagWell, EdgeBio).
7. Nylon cell strainers with 70  $\mu$ m pores.
8. Long-stem analytical funnel.
9. Pipet-Aid.
10. Sterile 10 mL plastic serological pipettes.
11. Microfuge tubes, 1.5 mL (e.g., Eppendorf).
12. PCR tubes, 0.2 mL.
13. Falcon tubes, 15 mL and 50 mL.
14. Nutator platform rotator.
15. Hemocytometer (e.g., Hausser Bright Line hemocytometer, Fisher Scientific).
16. Microcentrifuge and refrigerated centrifuge with rotor for 15 mL tubes.
17. Cold room, 4 °C.
18. Molecular biology grade water.
19. Sterile disposable filter unit, 500 mL.
20. Sterile 0.2  $\mu$ m syringe filter.
21. Sterile 10 mL plastic syringe.
22. Thermal cycler.
23. Real-time PCR machine.

24. A 64-bit computer with at least 1 TB hard disk and 16 Gb of memory for ATAC-seq data analysis.
25. Fluorescent microscope equipped for GFP detection and a 40× magnification objective or more.

## **2.2 Stock Solutions and Reagents**

All solutions are done in sterile, ultrapure purified water unless indicated.

1. Complete, EDTA-free Protease Inhibitors (e.g., Roche).
2. Stock solution of 2 M spermidine. Prepare by dissolving 2.904 g spermidine powder in 10 mL water. Aliquot 1 mL per 1.5 mL Eppendorf tube and store at  $-20^{\circ}\text{C}$ .
3. Stock solution of 200 mM spermine. Prepare by dissolving 0.4047 g spermine powder in 10 mL of water. Aliquot 1 mL per 1.5 mL Eppendorf tube and store at  $-20^{\circ}\text{C}$ .
4. Stock solution (1 L) of incomplete Nuclei Purification Buffer (NPBi): 20 mM MOPS, 40 mM NaCl, 90 mM KCl, 2 mM EDTA, 0.5 mM EGTA, adjusted to pH 7 with 2 M KOH. Filter-sterilize the solution and degas under vacuum for 10 minutes. Store at  $4^{\circ}\text{C}$  for up to 3 months.
5. Stock solution of 10% Triton X-100.
6. Stock solution of 10× DAPI. Prepare by dissolving 10 mg DAPI powder in 5 mL water, for a final concentration of  $2\text{ }\mu\text{g}/\mu\text{L}$ . Filter-sterilize the solution and store at  $4^{\circ}\text{C}$  in the dark for several months. To stain nuclei with DAPI, dilute the 10× DAPI solution to 1× using water (final concentration of  $0.2\text{ }\mu\text{g}/\mu\text{L}$ ), and use within 2–3 h.

## **2.3 Purification of Tagged Nuclei Using INTACT**

1. Plant material: tissue from transgenic plants expressing both NTF and BirA in the cell type of interest. INTACT transgenic lines targeting the root epidermal hair and non-hair cell types, as well as INTACT plasmid vectors are available from the Arabidopsis Biological Resource Center at Ohio State University.
2. Magnetic beads coupled to Streptavidin, 2–3  $\mu\text{m}$  diameter (e.g., M-280 Streptavidin Dynabeads, Life Technologies).
3. Protease inhibitor cocktail (e.g., Complete protease inhibitor, Roche).
4. Nuclei Purification Buffer (NPB): 20 mM, 40 mM NaCl, 90 mM KCl, 2 mM EDTA, 0.5 mM EGTA, 0.5 mM spermidine, 0.2 mM spermine, and 1× Roche Complete protease inhibitors, adjusted to pH 7 with 2 M KOH. Prepare by adding spermidine, spermine, and the protease inhibitors to NPB just before starting the INTACT nuclei purification procedure. Keep solution on ice, and use within 1 h of preparation.

5. Nuclei Purification Buffer containing 0.1% Triton X-100 (NPBt): 20 mM MOPS pH 7, 40 mM NaCl, 90 mM KCl, 2 mM EDTA, 0.5 mM EGTA, 0.5 mM spermidine, 0.2 mM spermine, and 0.1% (v/v) Triton X-100. Prepare by adding spermidine, spermine, and Triton X-100 to NPBi just before starting the INTACT nuclei purification procedure. Keep solution on ice, and use within 1 day of preparation.

#### **2.4 Purification of Total Nuclei Using Sucrose Sedimentation**

1. Plant material: fresh or frozen plant tissue.
2. Stock solution of 1 M Tris-HCl pH 8.
3. Stock solution of 1 M MgCl<sub>2</sub>.
4. Stock solution of 2 M sucrose.
5. Nuclei Purification Buffer (NPB): 20 mM MOPS pH 7, 40 mM NaCl, 90 mM KCl, 2 mM EDTA, 0.5 mM EGTA, 0.5 mM spermidine, 0.2 mM spermine, and 1× Roche Complete protease inhibitors. Prepare by adding spermidine, spermine, and Roche Complete protease inhibitors to NPBi just before starting the nuclei purification procedure. Keep solution on ice, and use within 1 h of preparation.
6. Nuclei Extraction Buffer 2 (NEB 2): 0.25 M Sucrose, 10 mM Tris-HCl pH 8, 10 mM MgCl<sub>2</sub>, 1% Triton X-100, and 1× Roche Complete Protease Inhibitors. Prepare solution just before use, keep on ice, and use within 1 h of preparation.
7. Nuclei Extraction Buffer 3 (NEB 3): 1.7 M Sucrose, 10 mM Tris-HCl pH 8, 2 mM MgCl<sub>2</sub>, and 0.15% Triton X-100, 1× Roche Complete Protease Inhibitors. Prepare solution just before use, keep on ice, and use within 1 h of preparation.

#### **2.5 Tagmentation of Chromatin by Tn5 Transposase**

1. Tagmentation-based library preparation kit (e.g., Nextera, Illumina, or equivalent).
2. PCR purification kit (e.g., MiniElute, Qiagen).
3. Nuclease-free water (e.g., Sigma, Ambion).
4. Nucleic acid decontamination detergent (e.g., DNA AWAY, or equivalent).

#### **2.6 Sequencing Library Preparation**

1. ATAC Primer 1 (AATGATACGGCGACCACCGAGATCTACACTCGTCGGCAGCGTCAGATGTG).
2. ATAC barcoded Primer 2 (CAAGCAGAAGACGGCATACGAGATNNNNNNNNGTCTCGTGGGCTCGGAGATGT); N's indicate the 8-base index sequence. Each library to be pooled for sequencing should be amplified with a different barcoded primer 2. See Supplementary Table 1 of [7] for all primer sequences.
3. High-Fidelity PCR master mix for library amplification (e.g., NEBNext High-Fidelity 2× PCR Master Mix, NEB).

4. Nucleic acid fluorescent dyes for qPCR (e.g., Solution of 20× EvaGreen dye, Biotium).
5. Passive fluorescence dye for qPCR normalization (e.g., Solution of 50× ROX dye, Invitrogen).
6. PCR Purification kit (e.g., MinElute, Qiagen).
7. PCR purification magnetic beads (e.g., Agencourt AMPure XP, Beckman Coulter).
8. 100% ethanol.
9. Horizontal electrophoresis gel box and power source.
10. A 302 nm ultraviolet transilluminator.
11. Library Quantification kit (e.g., NEBNext, NEB).

---

### 3 Methods

Users should either begin at Subheading 3.1 for affinity purification of nuclei using INTACT, or at Subheading 3.2 for isolation of total nuclei. In either case, the purified nuclei are used for tagmentation by Tn5 transposase in Subheading 3.3. All procedures are carried out at room temperature (25 °C) unless otherwise specified.

#### **3.1 Purification of Tagged Nuclei Using INTACT**

Purification of nuclei by INTACT gives access to cell type-specific nuclei and provides the highest purity nuclei with the least organelle contamination. However, established transgenic lines are required.

1. Excise the tissue (3 g of roots or 0.5 g of leaves) into a weigh boat on ice until the desired amount has been collected and weighed.
2. Grind the tissue to a fine powder in liquid nitrogen using a mortar and pestle.
3. Using a nitrogen-cooled metal lab spoon, quickly transfer the frozen tissue powder to another mortar containing 10 mL of ice-cold Nuclei Purification Buffer (NPB).
4. Thoroughly resuspend the powder in NPB by grinding it with a new, clean pestle (*see Note 1*).
5. Place a 70 µm nylon cell strainer in the center of a long-stemmed funnel and above a 15 mL tube on ice.
6. Use a 10 mL serological pipette to draw up the tissue suspension and filter it through the strainer. Collect the flow-through in the 15 mL tube.
7. Spin down the nuclei at  $1,200 \times g$  for 10 min at 4 °C.
8. Use a 10 mL serological pipet and then a 1 mL pipette tip as needed to carefully remove as much of the supernatant as possible without disturbing the pellet.

9. Gently resuspend the pellet in 1 mL of ice-cold NPB. Transfer the crude nuclei suspension to a 1.5 mL tube. Keep on ice.
10. Streptavidin beads preparation: Wash the appropriate amount of magnetic beads (25  $\mu$ L for nuclei from 3 g of roots or 10  $\mu$ L for 0.5 g of leaves) with 1 mL of ice-cold NPB in a 1.5 mL tube. Collect the beads on the magnetic rack. Discard the supernatant and resuspend the beads with ice-cold NPB to their original volume (e.g., 25  $\mu$ L). Keep on ice.
11. Bind the biotinylated nuclei to the streptavidin beads: Add the freshly prepared beads to the 1 mL of resuspended nuclei from **step 9**. Rotate on a nutator in a 4 °C cold room for 30 minutes. Work in the 4 °C cold room for **steps 12–23**.
12. Transfer the 1 mL bead–nuclei mixture to a 15 mL tube and slowly add 13 mL of ice-cold NPBt. Mix gently and place on a nutator for 30 s.
13. Place the 15 mL tube in the magnetic rack for 15 mL tubes for 2 min to capture the bead-bound nuclei along the walls of the tube.
14. Slowly remove the NPBt supernatant with a serological pipette, making sure not to disturb the beads on the side walls of the tube.
15. Gently resuspend the beads with 14 mL of ice-cold NPBt, mix gently, and place on a nutator for 30 s.
16. Place the 15 mL tube in the magnetic rack for 2 min to capture the nuclei and beads.
17. Repeat **steps 14** through **16** one more time, for a total of three washes.
18. Slowly remove the NPBt supernatant with a serological pipette.
19. Resuspend the beads in 1 mL of ice-cold NPBt.
20. Remove 25  $\mu$ L of this nuclei–bead suspension and place in a 0.6 mL tube on ice for later quantification of captured nuclei with a hemocytometer (*see step 24*).
21. Transfer the remaining nuclei–bead suspension to an ice-cold 1.5 mL tube.
22. Place the 1.5 mL tube in the magnetic rack for microfuge tube to capture the beads along the walls of the tube.
23. Carefully remove the NPBt supernatant and resuspend the bead-bound nuclei in 20  $\mu$ L of ice-cold NPB. Keep on ice until the nuclei are counted and ready for tagmentation (*see Note 2*).
24. To view and quantify nuclei under a light microscope, add 1  $\mu$ L of diluted DAPI solution (0.2  $\mu$ g/ $\mu$ L) to each 25  $\mu$ L aliquot of nuclei from **step 20**. Mix well, and place on ice for 5 min in the dark.



25. Use a hemocytometer to count the DAPI-stained, bead-bound nuclei and determine the total yield. Purified nuclei should appear as shown in Fig. 1c (*see* **Note 3**).
26. Use the calculated total yield to determine the volume of resuspended nuclei from **step 23** needed to obtain 50,000 nuclei for the ATAC-seq reaction.
27. Transfer the volume containing 50,000 nuclei to a new 0.2 mL tube.
28. Place the 0.2 mL tube in the 96-well magnetic separator plate to capture the beads along the walls of the tube.
29. Remove the supernatant, and resuspend the bead-bound nuclei in 50  $\mu$ L of ice-cold transposition reaction mix (*see* Subheading 3.3).
30. *Immediately* proceed to Subheading 3.3.

### 3.2 Purification of Total Nuclei Using Sucrose Sedimentation

This method for purification will provide nuclei suitable for ATAC-seq from *Arabidopsis* leaf or root tissue, and may also be effective for other plant species with modifications. Users should be aware that substantial organelle contamination will exist in nuclei samples prepared with sucrose sedimentation. As such, the sequencing depth of the ATAC-seq libraries will need to be deeper than those from INTACT-purified samples in order to account for this.

1. Excise the tissue (0.1–1 g of plant tissue) into a weigh boat on ice until the desired amount has been collected and weighed.
2. Grind the plant tissue to a fine powder in liquid nitrogen using a mortar and pestle (*see* **Note 4**).
3. Using a nitrogen-cooled metal lab spoon, quickly transfer the frozen tissue powder to another mortar containing 10 mL ice-cold Nuclei Purification Buffer (NPB).
4. Thoroughly resuspend the powder in NPB by grinding it with a new, clean pestle.
5. Place a 70  $\mu$ m nylon cell strainer in the center of a long-stemmed funnel and above a 15 mL tube on ice.
6. Use a 10 mL serological pipette to draw up the tissue suspension and filter it through the strainer. Collect the flow-through in the 15 mL tube.
7. Spin down the nuclei at  $1,200 \times g$  for 10 min at 4 °C.
8. Use a 10 mL serological pipet and then a 1 mL pipette tip to carefully remove as much of the supernatant as possible without disturbing the pellet.
9. Gently resuspend the pellet in 1 mL of ice-cold Nuclei Extraction Buffer 2 (NEB2). Transfer this suspension to a new 1.5 mL tube.

10. Spin down the resuspended nuclei at  $12,000 \times g$  for 10 min at 4 °C.
11. Carefully remove the supernatant and resuspend the pellet thoroughly in 300  $\mu$ L of Nuclei Extraction Buffer 3 (NEB3).
12. Add 300  $\mu$ L of ice-cold NEB3 to a new 1.5 mL tube.
13. Carefully layer the resuspended pellet from **step 11** on top of the fresh NEB3 from **step 12**.
14. Spin down the two layers at  $16,000 \times g$  for 10 min at 4 °C (*see Note 5*).
15. Carefully remove the supernatant and resuspend the nuclei pellet in 1 mL of cold NPB. Keep on ice.
16. Remove 25  $\mu$ L of this nuclei suspension and move to a fresh 0.6 mL tube on ice for quantification of isolated nuclei with a hemocytometer (*see step 18*).
17. Transfer the remaining isolated nuclei to an ice-cold 1.5 mL tube.
18. To view and quantify nuclei under a light microscope, add 1  $\mu$ L diluted DAPI solution (0.2  $\mu$ g/ $\mu$ L) to each 25  $\mu$ L of nuclei from **step 16**. Mix well, and place on ice for 5 min in the dark.
19. Use a hemocytometer to quantify the DAPI-stained nuclei and determine the total yield. Purified nuclei should appear as shown in Fig. 1c (*see Note 6*).
20. Use the calculated total yield to determine the volume of resuspended nuclei from **step 17** needed to obtain 50,000 nuclei for the ATAC-seq reaction.
21. Transfer the volume containing 50,000 nuclei to a new 0.2 mL tube.
22. Spin down the nuclei at  $1,500 \times g$  for 7 min at 4 °C.
23. Remove the supernatant, and resuspend the nuclei in 50  $\mu$ L of ice-cold transposition reaction mix (*see Subheading 3.3*).
24. *Immediately* proceed to Subheading 3.3.

### 3.3 Tagmentation with Tn5 Transposase

This step will fragment the chromatin while adding the adapter sequences needed for high-throughput sequencing.

1. Prepare the transposition reaction master mix in a 0.2 mL PCR tube on ice according to Table 1 and mix well. The volumes given in Table 1 are for a single reaction with 50,000 nuclei.
2. Place the nuclei resuspended in transposition reaction mix – either from Subheading 3.1, **step 29** if the nuclei were isolated using magnetic beads or from Subheading 3.2, **step 23** if the nuclei were isolated using sucrose sedimentation – in a thermal cycler block prewarmed to 37 °C.

**Table 1**  
**Transposition reaction mix**

Component	Volume ( $\mu$ L)
2 $\times$ TD buffer	25
Water	22.5
TDE1 Transposase	2.5
Total	50

3. Incubate for 30 min, gently mixing the reaction by hand every 5 min.
4. Purify the transposed DNA using the PCR purification kit according to the manufacturer's instructions.
5. Elute DNA in 11  $\mu$ L of elution buffer (EB) provided in the kit. DNA can now be stored at  $-20^{\circ}\text{C}$  until future use, or used immediately for PCR amplification.

### 3.4 PCR Amplification of the DNA Library

After tagmentation, this step will increase the abundance of the library fragments while also allowing the libraries to be barcoded.

1. Prepare the PCR amplification mix in a 0.2 mL tube on ice according to Table 2.
2. Mix well, and perform PCR cycling as described in Table 3 (*see Note 7*).
3. Once the thermal cycler reaches  $4^{\circ}\text{C}$ , remove the samples and place them on ice.
4. To determine the number of additional PCR cycles needed to adequately amplify the DNA library, prepare the qPCR Library Amplification Mix described in Table 4 in a 0.2 mL PCR tube. Keep the mixture on ice.
5. Perform thermal cycling in the qPCR machine according to Table 5.
6. To determine the optimal number of cycles needed to amplify the remaining 45  $\mu$ L of each library from **step 2**, view the linear fluorescence versus cycle number plot on the qPCR machine once the reaction is finished. The cycle number at which the fluorescence for a given reaction is at  $1/3$  of its maximum is the number of additional cycles ( $N$ ) that each library requires for adequate amplification (*see Note 8*).
7. Run the remaining 45  $\mu$ L of each PCR reaction from **step 3** according to Table 6.

**Table 2**  
**Transposed DNA amplification mix**

Component	Volume ( $\mu$ L)
Transposed DNA (from Subheading 3.3, step 5)	10
Water	10
25 $\mu$ M ATAC Primer 1	2.5
25 $\mu$ M ATAC barcoded Primer 2 <sup>a</sup>	2.5
2 $\times$ NEBNext High Fidelity PCR Mix	25
Total	50

<sup>a</sup>A different barcoded Primer 2 should be used for each library that is to be pooled into a single sequencing run.

**Table 3**  
**Thermal cycling conditions for transposed DNA amplification**

Cycle number	Temperature ( $^{\circ}$ C)	Time
1	72	5 min
	98	30 s
5 cycles	98	10 s
	63	30 s
	72	1 min
	4	Hold

**Table 4**  
**qPCR library amplification mix**

Component	Volume ( $\mu$ L)
Amplified library (from Subheading 3.4, step 3)	5
Water	0.45
25 $\mu$ M ATAC Primer 1	0.5
25 $\mu$ M ATAC barcoded Primer 2	0.5
20 $\times$ Evagreen dye	0.75
50 $\times$ ROX dye <sup>a</sup>	0.30
2 $\times$ NEBNext High Fidelity PCR Mix	7.5
Total	15

<sup>a</sup>ROX concentration may vary depending on the qPCR instrument. The amount described here is optimized for the ABI Step-One-Plus instrument.


**Table 5**  
**qPCR cycling conditions to determine additional library amplification cycles**

Cycle number	Temperature (°C)	Time
1	98	30 s
20 cycles	98	10 s
	63	30 s
	72	1 min

**Table 6**  
**Final library amplification**

Cycle number	Temperature (°C)	Time
1	98	30 s
<i>N</i> cycles	98	10 s
	63	30 s
	72	1 min
	4	Hold

8. Purify the libraries by mixing PCR purification magnetic beads with the reaction products at a 1.5:1 volume ratio of beads: PCR sample (*see Note 9*).
9. Incubate at room temperature for 5 min.
10. Place the 0.2 mL tube on the 96-well magnetic separator plate for 1 min to capture the PCR purification magnetic beads, and discard the supernatant.
11. With the tubes still in the magnetic plate, wash the beads twice for 30 s each with 200  $\mu$ L of 80% ethanol, without disturbing the bead pellet.
12. After the last wash, allow the beads to dry for 5 min to remove all traces of ethanol (*see Note 10*).
13. Remove the tube from the magnet and resuspend the bead pellet in 20  $\mu$ L 10 mM Tris pH 8.
14. Incubate at room temperature for 2 min.
15. Capture the beads on the magnet, and transfer the supernatant into a fresh 0.2 mL PCR tube on ice.
16. A small aliquot of the library, 1–2  $\mu$ L, can be run on a 2% agarose gel to visualize the abundance and size distribution of amplified libraries (Fig. 2a) (*see Note 11*).



**Fig. 2** ATAC-seq library preparation and high-throughput sequencing. **(a)** An amplified ATAC-seq library purified with AMPure XP beads (*lane “1”*) was resolved in a 2% agarose gel stained with ethidium bromide. *Lane “M”* is the molecular weight marker lane. Amplified library fragments generally range in size from 180 bp to several kb in size. The size distribution of the resolved gel may vary somewhat, but the final product should be free of adapter dimers (distinct band around 125 bp) and primer dimers (distinct band around 80 bp) (*see Note 11*). **(b)** Insert sizes of ATAC-seq paired-end reads from 50,000 nuclei isolated by INTACT from non-hair cells calculated using the InsertSizeMetrics option from Picard Tools (*see Note 13*). The distribution shows periodicity of the helical pitch of DNA for fragments smaller than 200 bp. Fragments containing one or more nucleosomes, with insert sizes of 150 or 300 bp were not observed using the transposase: nuclei and bead: DNA ratios described in this protocol. **(c)** Integrated Genome Viewer snapshot of four different libraries sequenced on the Illumina platform. The tracks shown are of ATAC sequencing reads from INTACT isolated nuclei from root hair cells (*orange*), root non-hair cells (*purple*), root tip tissue (*cyan*), and sucrose sedimentation isolated nuclei from 1 cm root tip tissue (*navy*). Gene tracks are shown below the ATAC-seq tracks and a 25 kb scale bar is shown

17. The purified libraries can now be stored at  $-20^{\circ}\text{C}$ .
18. Quantify the molar concentrations of the libraries using the Library Quantification kit according to manufacturer's directions.
19. Once quantified, the libraries are ready for pooling and high-throughput sequencing on a next-generation sequencing platform (*see Note 12*).

20. Analysis of the quality of the sequencing reads, alignment to the genome, examination of the fragment size distribution (Fig. 2b), and further downstream analyses can be performed as described briefly in **Note 13**. A genome browser shot of the typical *Arabidopsis* ATAC-seq data from libraries made using the procedures described here can be seen in Fig. 2c.

---

## 4 Notes

1. This protocol is optimized for 3 g of root or 0.5 g of leaf tissue from *Arabidopsis*. Ground leaf tissue contains more debris relative to roots, and therefore necessitates a lower amount of starting material to obtain highly purified nuclei. INTACT may also be performed on fresh tissue by chopping the tissue in NPB as opposed to grinding to a fine powder using liquid nitrogen. However, this approach does require the use of fresh tissue. The number of samples that can be run through INTACT purification simultaneously is mainly limited by the capacity of the magnetic rack for 15 mL tubes used for nuclei capture. For example, up to four separate samples can be processed in parallel using one DynaMag 15 magnetic rack. Using an INTACT line with nuclei labeled in the root epidermal non-hair cell type, approximately 200,000 purified nuclei can be obtained from 3 g of roots. Larger amounts of tissue can be used for purifying nuclei from less abundant cell types, and this generally only requires adjustments to the amount of streptavidin beads used and the volume of solution used for bead capture. *See* [11] for more details on variations in the INTACT procedure.
2. After isolating the bead bound nuclei, keep the sample on ice while quantifying the nuclei from the aliquot in Subheading 3.1, **step 23**. *Do not* freeze the isolated nuclei before doing fragmentation and library preparation. Freezing and thawing of isolated nuclei can disrupt protein–DNA interactions.
3. After DAPI staining, nuclei purified by INTACT can be easily identified and counted using a hemocytometer. The ideal setup for visualizing nuclei is under a mix of dim white light and DAPI channel fluorescence. The dim white light allows for visualization of the hemocytometer grid and the beads, and the DAPI fluorescence allows for the visualization of nuclei. A sample image of isolated bead-bound nuclei is shown in Fig. 1c. A nucleus is identified as a punctate circle with strong DAPI fluorescence that has several beads clustered around it. Minimal cellular debris or contaminating unbound nuclei should be observed in the final product. These contaminants may be further reduced by using fewer beads and by increasing

the volumes of NPB and NPBt used during purification as described in **Note 1**. We have successfully used as few as 20,000 to as many as 200,000 INTACT-purified nuclei in this procedure without altering any other parameters of the protocol presented here.

4. This protocol is optimized for less than 1 g of root or 0.5 g of leaf tissue. Ground leaf tissue contains more debris relative to roots, and therefore necessitates a lower amount of starting material to obtain purified nuclei. As with the INTACT protocol, sucrose sedimentation of nuclei may also be performed on fresh tissue by chopping the tissue in NPB as opposed to grinding to a fine powder using liquid nitrogen. However, this approach does require the use of fresh tissue. We recommend starting with the minimum amount of tissue needed to obtain the required number of nuclei (e.g., 50,000 per ATAC-seq reaction).
5. Proper separation of nuclei from other cellular debris requires the nuclei to pass through the sucrose cushion during centrifugation. The NEB3 resuspended nuclei should therefore be placed gently on top of NEB3 layer present in the tube. After centrifugation, the contaminating organelles and debris may be visible at the top of the tube. If leaf tissue was used, the top layer will become greener after centrifugation and the pellet will become noticeably less green than it was prior to centrifugation.
6. After DAPI staining, nuclei purified by sucrose sedimentation can be identified and quantified using a hemocytometer. A mixture of DAPI-channel fluorescence and white light illumination allows the stained nuclei and the hemocytometer grid to be seen simultaneously. A sample image of isolated nuclei is shown in Fig. 1c. A nucleus is identified as a punctate circle with strong DAPI fluorescence. The nucleus is typically ~5  $\mu\text{m}$  in size and can be easily identified at 200 $\times$  and 400 $\times$  magnifications. Cellular debris may be observed in the final preparation, but this generally does not affect the outcome of the ATAC-seq procedure. To reduce cellular debris contamination, starting tissue can be chopped with a razor blade (*see Note 4*) and/or additional NEB3 wash steps may also be done by repeating Subheading 3.2, steps 11–14 for a second sucrose cushion centrifugation.
7. Ensure that all work surfaces, pipettes, and reagents needed for amplification and library preparation are free of DNA contamination by wiping them down with a nucleic acid decontamination detergent 10 min before starting work. For library amplification, unique barcoded primers are used for each sample if multiple libraries are to be sequenced in an individual flow cell lane. The sequences of all primers can be found in the supplementary material of [7].



8. The number of PCR cycles needed to amplify ATAC libraries is determined by the PCR reaction in Subheading 3.4, step 6. We recommend using the minimum number of cycles necessary to obtain a sufficient molar amount of library for Illumina sequencing. This must be determined empirically and will also depend on the number of libraries to be pooled for sequencing. Typically, when the libraries are amplified to 1/3 of their maximum fluorescence detected by qPCR the quantified molarities are from 50 to 300 nM. The variation in the final product is dependent on how reliably nuclei were counted in Subheading 3.1, step 25 or Subheading 3.2, step 19.
9. The ratio of the PCR purification magnetic beads to the amplified library volume determines the size of the purified DNA fragments that are isolated. The 1.5 bead-to-amplified library reaction ratio results in the isolation of DNA fragments shown in Fig. 2a. Using ratios that have higher proportions of beads may result in purification of sequencing adapters and PCR primers, which can negatively affect sequencing.
10. A drying time of 5 min is generally sufficient to remove all traces of ethanol from the beads, but this time may vary based on humidity and room temperature. Ensure that all ethanol has evaporated before moving on to the next step. Do not allow beads to dry to the extent that the pellet begins to crack as this will decrease the purification yield.
11. Libraries can generally be visualized by agarose gel electrophoresis followed by ethidium bromide staining (Fig. 2b). Sensitivity can be greatly increased by staining the gel with a nucleic-acid specific green fluorescent dye (e.g., SYBR green stain, Qiagen, or equivalent) or using microfluidics-electrophoresis trace analysis (e.g., Agilent Bioanalyzer or equivalent instrument), if available. The libraries that we have prepared using this method generally present as a DNA smear starting at ~180 bp and ranging to greater than 1 kb, with peak intensity falling between ~180 and 500 bp (*see* Fig. 2a). The original publication on ATAC-seq [7] reported a nucleosome-like periodicity in the library size distribution, but we have not observed this phenomenon either by electrophoresis or by the estimation of fragment size distribution based on the distance between paired-end sequencing reads, as shown in Fig. 2b. This lack of observed nucleosome fractions may be due to the size selection of library fragments by PCR purification magnetic beads and the low transposase-to-nuclei ratio described in this protocol.
12. Paired-end sequencing is recommended in order to maximize the number of transposase integration events that can be observed in a given sample, and to allow measurement of the

length of the sequenced fragments (Fig. 2b). To identify open chromatin regions in *Arabidopsis*, users should aim to obtain at least 10–20 million reads per library that map to the nuclear genome. For transcription factor footprinting the number of nuclear genome-mapping reads should be increased to at least 60 million per library [12]. When using sucrose sedimentation for nuclei purification, users should expect ~50% of reads to map to the nuclear genome, while the use of INTACT purification will increase this number to >90%.

13. Sequencing reads are checked for overall quality using FastQC (<http://www.bioinformatics.babraham.ac.uk/projects/fastqc/>) or equivalent. The reads are aligned to the TAIR10 *Arabidopsis thaliana* genome ([https://www.arabidopsis.org/download/index-auto.jsp?dir=%2Fdownload\\_files%2FGenes%2FTAIR10\\_genome\\_release](https://www.arabidopsis.org/download/index-auto.jsp?dir=%2Fdownload_files%2FGenes%2FTAIR10_genome_release)) or the most recent genome annotation version using Bowtie2 (<http://bowtie-bio.sourceforge.net/bowtie2/index.shtml>). The resulting SAM file is converted to a binary BAM file, which is sorted and indexed using Samtools (<http://samtools.sourceforge.net/>). The quality of the resulting BAM file, including fragment size distribution, is analyzed using Picard Tools (<https://broadinstitute.github.io/picard/>). Alignment data is visualized using the Integrated Genome Viewer (<http://software.broadinstitute.org/software/igv/>). For ease of visualization, BAM files were converted to BigWig files using DeepTools BamPECoverage tool (<http://deeptools.readthedocs.io/en/latest/index.html>). Further downstream analyses of ATAC-seq data include calling peaks with HOMER (<http://homer.salk.edu/homer/index.html>), editing BED files with bedtools (<http://bedtools.readthedocs.io/en/latest/>), and identifying transcription factor footprints using pyDNase (<http://pythonhosted.org/pyDNase/>).

---

## Acknowledgments

This work was supported by the National Science Foundation Grant No. IOS-1238243. We thank Paja Sijacic and Shannon Torres for helping to optimize the protocol for nuclei isolation and for suggestions on the manuscript.

## References

1. Li B, Carey M, Workman JL (2007) The role of chromatin during transcription. *Cell* 128:707–719
2. Pajoro A, Muiño JM, Angenent GC, Kaufmann K (2017) Profiling nucleosome occupancy by MNase-seq: experimental protocol and computational analysis. In: Bemer M, Baroux C (eds) *Plant chromatin dynamics: methods and protocols*. Springer, New York, NY. doi:10.1007/978-1-4939-7318-7\_11

3. Song L, Crawford GE (2010) DNase-seq: a high-resolution technique for mapping active gene regulatory elements across the genome from mammalian cells. *Cold Spring Harb Protoc* 2010:pdb.prot5384. doi:[10.1101/pdb.prot5384](https://doi.org/10.1101/pdb.prot5384)
4. Ken Z (2005) Micrococcal nuclease analysis of chromatin structure. *Curr Protoc Mol Biol* 21 (1):1–17
5. Desvoyes B, Vergara Z, Sequeira-Mendes JO, Madeira S, Gutierrez C (2017) A rapid and efficient ChIP protocol to profile chromatin binding proteins and epigenetic modifications in bulk *Arabidopsis* tissue. In: Bemer M, Baroux C (eds) *Plant chromatin dynamics: methods and protocols*. Springer, New York, NY. doi:[10.1007/978-1-4939-7318-7\\_5](https://doi.org/10.1007/978-1-4939-7318-7_5)
6. Park PJ (2009) ChIP-seq: advantages and challenges of a maturing technology. *Nat Rev Genet* 10:669–680
7. Buenrostro JD, Giresi PG, Zaba LC, Chang HY, Greenleaf WJ (2013) Transposition of native chromatin for fast and sensitive epigenomic profiling of open chromatin, DNA-binding proteins and nucleosome position. *Nat Methods* 10:1213–1218
8. Buenrostro JD, Wu B, Chang HY, Greenleaf WJ (2015) ATAC-seq: a method for assaying chromatin accessibility genome-wide. *Curr Protoc Mol Biol* 109:21.29.1–21.29.9
9. Gendrel A, Lippman Z, Martienssen R, Colot V (2005) Profiling histone modification patterns in plants using genomic tiling microarrays. *Nat Methods* 2:213–218
10. Deal RB, Henikoff S (2010) A simple method for gene expression and chromatin profiling of individual cell types within a tissue. *Dev Cell* 18:1030–1040
11. Wang D, Deal RB (2015) Epigenome profiling of specific plant cell types using a StreamLined INTACT protocol and ChIP-seq. *Methods Mol Biol* 1284:3–25
12. Lu Z, Hofmeister BT, Vollmers C, DuBois RM, Schmitz RJ (2016) Combining ATAC-seq with nuclei sorting for discovery of cis-regulatory regions in plant genomes. *Nucleic Acids Res*. doi:[10.1093/nar/gkw1179](https://doi.org/10.1093/nar/gkw1179)